

DIPHOTERINE[®] solution

CAPTURES THE CORROSIVE AND ELIMINATES IT



PREVOR

ANTICIPATE AND SAVE

Toxicology Laboratory & Chemical Risk Management

In accordance with the Standard EN 15154

CE 0459

The chemical burn

Two categories of chemicals can provoke a burn by a direct contact: corrosives and irritants.

The seriousness of the chemical burn depends on the number of cells destroyed and will be more significant in the case of corrosives. In some cases, the aggressive product can also have toxic or harmful effects, the consequences of which should not be neglected.

■ The aggressive chemicals

CORROSIVE



Concentrated acids and bases



SIGNIFICANT REACTIONS



IRREVERSIBLE EFFECTS

IRRITANT



Solvents, oils...



WEAK REACTIONS



REVERSIBLE EFFECTS



Specific danger if combined with toxic risk

(example - Hydrofluoric acid)


■ The mechanism of a chemical burn

THE EXCHANGE

SKIN and EYE


CHEMICAL

BURN



A chemical burn is caused by the capability of corrosives and irritants to create an exchange (ion, proton, electron...) with the tissues of the skin or the eye. The degree of the burn will depend on the number of molecules destroyed and on the type of modification (reversible or irreversible)

■ The action of soda NaOH:



Healthy cells as seen under a microscope

Addition of soda 0.5N (2%)

Cells are completely destroyed




The principles of emergency washing



To stop the aggressive reactions

- There are 6 types of aggressive chemical reactions




To be effective regardless of the type of accident and to avoid the risk of error for the victim, it is necessary to be able to stop these 6 aggressive reactions.



A POLYVALENT PRODUCT IS ESSENTIAL



To stop the evolution of the burn



A chemical burn is initiated by the contact between the aggressive product and the skin or the eye. Following this contact, a part of the aggressive chemical will penetrate into the tissue and cause destruction of the cells.

To effectively decontaminate the splash, it will be necessary, not only to decontaminate the product on the surface, but also to control its penetration inside the tissue.



A PRODUCT ABLE TO STOP THE AGGRESSIVE CHEMICAL IS ESSENTIAL

■ Factors influencing the penetration

- The length of exposure time
- The temperature
- The type of product and its concentration




First-aid treatment: From water...

■ The principles of washing with water


- a washing of the surface to remove quickly the aggressive product
- the dilution of the chemical to reduce its aggressiveness
- a universal product avoids the risks of error at the time of the accident

■ What are its limitations?


- the concentrated products which penetrate very quickly
- the intervention time of 10 seconds which is not always achievable
- the washing comfort : risk of hypothermia under a water shower
difficulty to open the eye



Healthy cells as seen under a microscope.



The beginning of washing: water penetrates inside the cells and makes them expand.




The end of washing : cells are destroyed because water causes them to burst.

*photo : ACTO, Pr. Norbert Schrage, Aix la Chapelle

Ex vivo test on the eye

Efficacy of the active washing with DIPHOTERINE® solution:
Ex Vivo EVEIT Model - OCT



Rabbit corneas, 16 minutes after an application of 500 µL 1M NaOH for 20 s.

- a) without any washing
- b) with DIPHOTERINE® solution washing

Corrosive penetration is stopped. The structural changes in the stroma are negligible. The endothelium is completely preserved.

Sources : Burgher, Mathieu, Fosse, Rihawi, Gérard, Merle, Schrage, Ocular chemical burn : Experimental proof of the influence of key parameters on the diffusion and the decontamination. 114th Congress of the SFO Paris May 2008.



The ideal product in case

Retain the advantages of water:

- Fast washing of the surface
- Single protocol

EXPERIENCE FEEDBACKS:

Result on the eye with a delayed washing



DIPHOTERINE® solution a better healing

> Ammonia ocular burn, washed

- Assessment before washing scale, which usually requires
- After a washing with 1 litre and 6 months of appropriate requiring a corneal graft, a

Sources : Gérard M, Merle H, Chiambaretta F, Rihawi, Ocular burn. Burns 2002;7:670-3.


DIPHOTERINE® solution is recognized as a reference solution by the SFO, French Society of Ophthalmology




...to DIPHOTERINE[®] solution

■ The principles of washing with DIPHOTERINE[®] solution


- It is a liquid which enables it to obtain the same effect as water on the surface of the skin or eye.
- DIPHOTERINE[®] solution is an amphoteric chelating agent, which enables it to stop the aggressiveness of the chemicals in a polyvalent way (*for HF and its derivatives, use the HEXAFLUORINE[®] solution*).
- DIPHOTERINE[®] solution mechanism can be illustrated as follows:



BASE ACID
DIPHOTERINE[®] solution will attract the chemical in contact with the tissues.




BASE ACID
The acid site of DIPHOTERINE[®] solution will fix bases to make them harmless.




BASE ACID
The basic site of DIPHOTERINE[®] solution will fix acids to make them harmless.


- DIPHOTERINE[®] solution, unlike water, will **stop the penetration** of the chemical.



Healthy cells as seen under a microscope.



Beginning of the washing with DIPHOTERINE[®] solution: the cells contract slightly.



End of the washing: the cells are preserved.

of a chemical splash must:

Bring improvements to compensate for the limitations of water:

- Guarantee a total efficacy whatever the product
- Increase the intervention time
- Improve washing comfort to increase the effectiveness

retains an interest in delayed washing for

ed with DIPHOTERINE[®] solution after 1 hour

with DIPHOTERINE[®] solution: grade IV burn on the Roper-Hall a corneal graft to achieve healing.

of DIPHOTERINE[®] solution within one hour after the splash e care, the victim regained a visual acuity of 14/20 without first in the history of medicine.

gal Rigal D, Schrage N. An amphoteric rinse used in the emergency treatment of a serious

RESULTS

EXPERIENCE FEEDBACK:

Results on the skin in emergency

Independent retrospective study conducted by Dr Donaghue, chief medical officer of Alcoa Australia (alumina refineries). It covers 180 cases studied from May 1st 2005 to April 30th 2008, specifically on the skin.

Severity scale and associated signs	FIRST AID EMERGENCY SOLUTION	
	First washing with a DIPHOTERINE [®] DAP	First washing with water
1 [no sign]	52,9% (73 cases)	21,4% (9 cases)
2 [erythema]	39,1% (54 cases)	54,8% (23 cases)
3 [blisters]	7,2% (10 cases)	19,0% (8 cases)
4 [more serious]	0,7% (1 case)	4,8% (2 cases)
TOTAL	100% (138 cases)	100% (42 cases)

Criteria used:

- 1 : Time elapsed between the chemical splash and clinical evaluation.
- 2 : Time elapsed between the chemical splash and application of DIPHOTERINE[®] solution.
- 3 : Percentage of skin surface affected by the chemical.



How to use

A sterile washing solution

Dispensers for the skin, in the form of micronised sprays (optimising the surface of contact)

ALL THESE PACKAGES MEET THE NEW EUR



MINI DAP



MICRO DAP



> To decontaminate an entire body

Typical equipment for locations where there are large quantities of chemicals :

- Production areas
- Warehousing areas
- Off-loading and decanting areas

> To decontaminate a face or an arm

An ideal complement for a Wall Mounted Eyewash, where the quantities of chemicals are limited:
Laboratories

> To decontaminate a hand

Ideal personal equipment for small splashes:
Workshops
Maintenance personnel

Washing must begin within the first minute following the splash



Protocol
DIPHOTERINE®



1/ Go away from the danger



2/ Get undressed

DIPHOTERINE[®] solution?

To be used immediately after the accident

Dispensers for the eyes, equipped with an ergonomic eyecup, which assists the opening of the eye for a more effective washing

EUROPEAN STANDARD EN 15154 PARTS 3 AND 4



> Kit for chemical splashes

In areas where fixed equipments are appropriate:

- Laboratories
- Production areas
- Warehouse areas
- Decanting areas

> To decontaminate an eye

To carry on oneself or to put in the work area:

- Laboratories
- Warehousing areas

Portable eyewash: when the equipment must be mobile: first aid kits, emergency vehicles, infirmaries

> Personal equipment to decontaminate an eye:

Ideal for the maintenance personnel
Imperatively worn by individual

Washing must begin within the first 10 seconds following the splash

Washing must begin within the first minute following the splash

For the first aid intervention with solution in the event of a chemical splash



3/ Wash as quickly as possible, respecting the instructions for use of DIPHOTERINE[®] solution



4/ Alert



5/ Seek medical advice





PREVOR

ANTICIPATE AND SAVE

Toxicology Laboratory & Chemical Risk Management

People saving People

**WE HELP YOU TO FACE
CHEMICAL HAZARDS**



References They are trusting us

AIR LIQUIDE - AKZO NOBEL - AOSTE - AP-HP - ARCELORMITTAL - ARJOWIGGINS
ARKEMA - ASTRAZENECA - BAYER - BONDUELLE - BRENNTAG - CARGILL - CEA
CEPHALON - CLARIANT - COCA COLA - CRISTAL UNION - DANONE - DASSAULT
DCNS - DOW - ECOLAB - EDF - ERAMET - FAREVA - GDF SUEZ - HOLCIM
HUNTSMAN - IBERDROLA - IMERYS - INTERNATIONAL PAPER - IPSEN - IRSN

ITALCEMENTI - LACTALIS - L'OREAL - NESTLE - PREZIOSO - QUARON - RIO TINTO
ROQUETTE - ROULLIER - SAFRAN - SAINT GOBAIN - SAIPOL - SANOFI - SERVIER
SIAAP - SIGMA ALDRICH - SNCF - SOLVAY - SUDZUCKER - TEREOS - TERRENA
THALES - TOTAL - VALEO - VALLOUREC ET MANNESMAN - VEOLIA - VINCI - YARA
YVES ROCHER...



PREVOR
ANTICIPATE AND SAVE

www.prevor.com

Moulin de Verville
95760 Valmondois
Tel. 00 33 (0)1 30 34 76 76
Fax. 00 33 (0)1 30 34 76 70