



Protocol for Hexafluorine® Solution Use

In cases of chemical splashes/contamination of fluorides within an acidic medium


ACTION WITHIN THE FIRST MINUTE



1 eye
Apply within 60 secs



Single Use
1 LPMF
500ml



Small area-Full Body
Apply within 60 secs



Use
1 DAPF
5L

Depending on the company's recommended medical protocol, apply a specific antidote such as calcium gluconate gel to the contaminated area.

DO NOT APPLY CALCIUM GLUCONATE TO THE EYES

In the case of delayed washing, the use of calcium gluconate gel or similar antidotes is particularly important. This is due to risk of systemic toxicity in these cases.

Please always refer Calcium Gluconate Gel Manufacturer's Instructions for Use.

In all cases, washing must be followed by medical examination.

Hexafluorine® Solution has limited efficacy on chemicals that are not acids.
Diphoterine® Solution is better adapted for this use.

HEXAFLUORINE® SOLUTION – WASHING PROTOCOL

1. Start washing uncovered areas within 60 seconds as per protocols overleaf.
2. Remove contaminated clothing and, if applicable, contact lenses.
3. Using full contents, continue washing all contaminated areas.
4. Do not interrupt washing process.
5. Do not put contaminated clothing back on.
6. Seek medical advice.

HEXAFLUORINE® SOLUTION – GENERAL INSTRUCTIONS

Never delay washing

- For optimal efficacy use Hexafluorine® Solution as the primary action.
- Use the entire contents of the container.
- If there is no Hexafluorine® Solution available use water and then wash with Hexafluorine® Solution as soon as possible.
- For a chemical contact time greater than 1 minute prolong the washing of the contaminated are with Hexafluorine® Solution for 3 to 5 times the duration of contact.
- In the case of a delayed ocular splash it is not necessary to continue washing with Hexafluorine® Solution for more than 15 minutes.

Always seek medical advice after decontamination

- After an active washing of the eye with Hexafluorine® Solution the use of **Afterwash II®** is recommended (after 15 minutes) to facilitate a more rapid return of the eye to a physiological state.
- If oral mucosa is affected by the splash, rinse the mouth with Hexafluorine® Solution, then spit out.
- If the ear canal is affected, wash as rapidly as possible with Hexafluorine® Solution by applying 500ml inside the canal, leaning the head to one side to allow the Solution to flow out of the ear

As in any case of unilateral washing of one ear with a liquid at room temperature, a dizzy feeling, without any serious consequences can occur, this will disappear after a few minutes.

Container	Average Diffusion Time
LPMF(500ml)	3 minutes
DAPF (5L)	5 minutes

Hexafluorine® Solution has limited efficacy on chemicals that are not acids.
Diphoterine® Solution is better adapted for this use.



DIPHEX SOLUTIONS LIMITED
Safety through Science
01622 851000
www.diphex.com